



3-Ply Disposable Masks

FDA Registered Manufacturer

Technical Data Sheet



DETAILS:

- Flat, easy on/easy off mask with standard pleats
- Easy to breathe through
- Bi-directional pleats for increased facial comfort and fit
- Flat, easy on/off ear loop mask with standard pleats.
- Lightweight, comfortable and easy to breathe through.
- To be used in environments with minimal fluid exposure.

General patient care and other environments with minimum fluid exposure, isolation, IC practices, Central processing, SPD, clinics, long term care facilities and home health care.

Face Masks and Respirators are Personal Protective Equipment (PPE) that can be used to protect against airborne droplets containing viruses and bacteria, such as those generated by coughs and sneezes.

Using Personal Protective Equipment (PPE's) can protect against infection from many viruses including Coronavirus (e.g. SARS, MERS, COVID-19) and Flu Viruses (e.g. H1NI, Swine Flu, Bird Flu)

Features

- FDA registered manufacturer
- Manufactured to meet stringent guidelines
- Flat, easy on/off ear loop mask
- Standard pleats
- Lightweight, comfortable and easy to breathe
- Intended for use in environments with minimal fluid exposure.

Materials

The following materials are used in the production of 3-ply masks:

Straps	Braided Polyisoprene
First Layer	Polypropylene spun-bonded non-woven fabric
Second Layer	Meltblown non-woven fabric
Third Layer	Polypropylene spun-bonded non-woven fabric

This mask does not contain components made from natural rubber latex.

Maximum mass of products = 11.3g

50 pcs



Use For

- General patient care and other environments with minimum fluid exposure, isolation, IC practices, Central processing, SPD, clinics, long term care facilities and home health care.
- Always follow *User Instructions* and use in manner as indicated

Do Not Use For

- **DO NOT** use for gases or vapors (i.e. anesthetic gases such as isoflurane or vapors from sterilants such as glutaraldehyde.)
- **DO NOT** use in any manner not indicated in the User Instructions

Approvals and Standards

- Surgical masks are cleared by the FDA for use as a surgical mask.
- Do not have the word NIOSH and the approval type (i.e. N95) printed on the product.
- Are not designed to protect your lungs from airborne hazards.
- Are typically donned for a specific procedure and disposed of afterward.
- Breathability <2.1
- Filtration BFE >99%
- Filtration PFE >98%

Packaging	Product Code	Order Code
50 Masks/box, 10 boxes/case	3Ply-MSK	ES3PMSK50

Fitting Instructions

**Must be followed each time the mask is worn.
Before fitting device, ensure hands are clean.**


All mask components should be inspected for damage prior to each use.

1. Pull the bottom of the face mask over your chin.
2. Hold the mask by the ear loops and place a loop around each ear.
3. Check mask tightness before going around others:
 - Cover your face with both hands and exhale quickly.
 - To remove the mask: Remove it from behind (do not touch the front of the mask); discard immediately in an appropriate disposal bin. Thoroughly clean hands with alcohol-based hand rub or soap and water.

Three layers of protection

Polypropylene spunbonded non-woven fabric
Effectively block visible objects such as droplets

Meltblown nonwoven
Filter non-oily particulate matter from the air



Polypropylene spunbonded non-woven fabric
Absorb the hot air exhaled from the body, make the skin dry and comfortable

Warnings and Use Limitations

Always be sure that the complete product is:

- Suitable for the application;
 - Fitted correctly;
 - Worn during all periods of exposure;
 - Replaced when necessary.
- It is recommended that fit testing be conducted before assigning a mask to an individual. If you cannot achieve a proper fit then do not enter contaminated area or be around others.
 - Inspect mask before each use to ensure it is in good working condition. Examine all the mask parts for signs of damage including the two straps and fabric. The mask should be disposed of immediately upon observation of damaged or missing parts.
 - Leave the contaminated area immediately and contact supervisor if dizziness, irritation or other distress occurs.
 - Dispose of used product in accordance with applicable regulations.
 - All masks should be used in accordance with local regulations.
 - Do not alter, repair, wash, and abuse or misuse the mask.
 - The mask can help protect the wearer's lungs against certain airborne contaminants; however, it will not prevent entry through other routes such as the skin or eyes, which would require additional personal protective equipment (PPE).
 - The mask is designed for occupational/professional use by adults who are properly trained in its use and limitations. The mask is not designed to be used by children.
 - Individuals with a compromised respiratory system, such as asthma or emphysema, should consult a physician and complete a medical evaluation prior to use.
 - Maximum Operating Temperature: +50 degrees Celsius.
 - The filtration efficiency of the mask may decrease in the presence of oily mists.

Storage and Transportation

Shelf life of the unopened product is two (2) years from date of manufacture when stored within temperature range of -20°C to +30°C and at less than 80% relative humidity. Before initial use, always check the product for signs of deterioration. When storing or transporting this product use original packaging provided.

Warning

Mask must not be used until you have determined whether usage will be in accordance with manufacturer's instructions. The wearer must be trained in the proper fitting and use of this product. Failure to follow all instructions and warnings on the use of this product and/or failure to wear this mask during all times of exposure can reduce mask effectiveness and result in illness or death. All masks should be used in accordance with United States standards or the country of origin. It is recommended that fit testing be conducted before assigning a mask to an individual. If you cannot achieve a proper fit, do not enter contaminated areas.

Important Notice

To the extent permitted by law, EnviroShield Supplies LLC ("EnviroShield") shall not be liable for any loss or damage including any loss of business, loss of profits, or for any indirect, special, incidental or consequential loss or damage arising from reliance upon any information herein provided by EnviroShield. Nothing in this statement will be deemed to exclude or restrict EnviroShield's liability for death or personal injury arising from its negligence.

Pricing

Pricing does not include shipping costs or any applicable taxes. Both ground shipping and air freight are available. Payment must be made in full prior to shipping. Escrow services are available but will incur an additional charge that must be paid by buyer.

Due to the volatility and high demand for masks, prices are subject to change without notice until payment is made in full.

Quantity	Price/Mask
500	\$1.00
1000	\$0.95
5000	\$0.90
10k	\$0.85
25k	\$0.80
50k	\$0.75
100k	\$0.70
250k	\$0.65
500k	\$0.63
1mm	\$0.60



Enviro Shield LLC dba PPEsite.org
16990 Dallas Pkwy Ste 200
Dallas, TX 75248

Helpline and Customer Service
(866) 568-1411

Email: orders@PPEsite.org
Web: PPEsite.org